

Anti-Hev b3 antibody [10] (STJA0002871)

STJA0002871

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Small rubber particle protein is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA, WB
Host/Source	Mouse
Reactivity	HEVB3

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	10
Concentration	
Conjugation	Unconjugated
Purification	Protein G purification
Dilution Range	ELISA 12000. Reacts as detection antibody in capture ELISA with Hev b3 binding monoclonal antibody A3-704-100. WB 12000 Monoclonal antibody working titer has to be established practically for each particular antigen and assay format
Formulation	PBS, 0.1% sodium azide
Isotype	IgG1
Storage Instruction	Store at -20 to -70 °C upon receipt. Divide antibody into aliquots prior usage. Avoid multiple freeze-thaw cycles as product degradation may result

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Hev b3-MBP fusion protein expressed and purified from E. coli. Hev b3
Immunogen Region	
Specificity	Recombinant and native Hev b3
Immunogen Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081