

Anti-Hev b5 antibody [4] (STJA0002872)

STJA0002872

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Primary antibodies |
| Short Description | Mouse monoclonal antibody anti-Major latex allergen is suitable for use in ELISA research applications. |
| Applications | ELISA |
| Host/Source | Mouse |
| Reactivity | HEVB5 |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | Monoclonal |
| Clone ID | 4 |
| Concentration | |
| Conjugation | Unconjugated |
| Purification | Protein G purification |
| Dilution Range | ELISA 132 000 Reacts as binding antibody in capture ELISA with Hev b5 detection monoclonal antibody A3-708-100. |
| Formulation | PBS, 0.1% sodium azide |
| Isotype | IgG1 |
| Storage Instruction | Store at -20 to -70 °C upon receipt. Divide antibody into aliquots prior usage. Avoid multiple freeze-thaw cycles as product degradation may result |

TARGET INFORMATION

| | |
|---------------------------|---|
| Gene ID | |
| Gene Symbol | |
| Uniprot ID | |
| Immunogen | Hev b5-MBP fusion protein expressed and purified from E. coli. Hev b5 |
| Immunogen Region | |
| Specificity | Recombinant and native Hev b5 |
| Immunogen Sequence | |