

Anti-Hev b6.02 antibody [2] (STJA0002874)

STJA0002874

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Hevein is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA, WB
Host/Source	Mouse
Reactivity	HEVB6, 02

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	2
Concentration	
Conjugation	Unconjugated
Purification	Protein G purification
Dilution Range	ELISA 164 000, WB. Reacts as binding antibody in capture ELISA with Hev b6.02 detection monoclonal antibody A3-710-100.
Formulation	PBS, 0.1% sodium azide
Isotype	IgG1
Storage Instruction	Store at -20 to -70 °C upon receipt. Divide antibody into aliquots prior usage. Avoid multiple freeze-thaw cycles as product degradation may result

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Hev b6.02-MBP fusion protein expressed and purified from E. coli. Hev b6.02
Immunogen Region	
Specificity	Recombinant and native Hev b6.02
Immunogen Sequence	