

Anti-GFRa-2 antibody (STJA0002878)

STJA0002878

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Chicken polyclonal antibody anti-GFRa-2 is suitable for use in ELISA, Western Blot and Immunofluorescence research applications.
Applications	ELISA, WB, IF
Host/Source	Chicken
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	GFRa-2 affinity purified
Dilution Range	WB 0.2 to 1 µg/ml; ELISA 6.25 to 8 ng/ml; IF 1 to 2 µg/ml Polyclonal antibody working titer has to be established practically for each particular antigen and assay format
Formulation	Concentrated ammonium sulphate in PBS pH 7.4
Isotype	IgY
Storage Instruction	Store at +4°C upon receipt. As product is (NH ₄) ₂ SO ₄ precipitate, mix well by pipetting or vortexing prior use

TARGET INFORMATION

Gene ID	2675
Gene Symbol	GFRA2
Uniprot ID	GFRA2_HUMAN
Immunogen	Recombinant His-tagged human GFR Alpha-2 protein produced using CHO-based proprietary suspension cell line. For production of hGFR Alpha-2, glycosylphosphatidylinositol GPI-anchor was removed and protein was secreted to the cell culture supernatant.
Immunogen Region	
Specificity	Human GFRa-2, doesn't cross-react with hGFRa-1, hGFRa-3, hGFRa-4
Immunogen Sequence	