

Anti-GFRa-3 antibody (STJA0002879)

STJA0002879

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Chicken polyclonal antibody anti-GFRa-3 is suitable for use in ELISA, Western Blot and Immunofluorescence research applications.
Applications	ELISA, WB, IF
Host/Source	Chicken
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	GFR Alpha-3 affinity purified
Dilution Range	WB 1 to 2 µg/ml; ELISA 3 to 6.25 ng/ml; IF 0.3 to 0.67 µg/ml. Polyclonal antibody working titer has to be established practically for each particular antigen and assay format
Formulation	Concentrated ammonium sulphate in PBS pH 7.4
Isotype	IgY
Storage Instruction	Store at +4°C upon receipt. As product is (NH ₄) ₂ SO ₄ precipitate, mix well by pipetting or vortexing prior use

TARGET INFORMATION

Gene ID	2676
Gene Symbol	GFRA3
Uniprot ID	GFRA3_HUMAN
Immunogen	Recombinant His-tagged human GFRA-3 protein produced using CHO-based proprietary suspension cell line. For production of hGFRA-3, glycosylphosphatidyl-inositol GPI-anchor was removed and protein was secreted to the cell culture supernatant. Protein w
Immunogen Region	
Specificity	Human GFRA-3, doesn't cross-react with hGFRA-1, hGFRA-2, hGFRA-4
Immunogen Sequence	