

Anti-Neurturin antibody (STJA0002881)

STJA0002881

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Neurturin is suitable for use in ELISA, Western Blot and Immunofluorescence research applications.
Applications	ELISA, WB, IF
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Immunogen affinity purification
Dilution Range	WB 0.2 to 1 µg/ml; ELISA 0.01 to 0.02 µg/ml; IF 0.5 to 2 µg/ml Antibody doesn't cross-react to other GDNF-family proteins Polyclonal antibody working titer has to be established practically for each particular antigen and assay format
Formulation	Concentrated ammonium sulphate in PBS pH 7.4
Isotype	IgG
Storage Instruction	Store at +4°C upon receipt. As product is (NH ₄) ₂ SO ₄ precipitate, mix well by pipetting or vortexing prior use

TARGET INFORMATION

Gene ID	4902
Gene Symbol	NRTN
Uniprot ID	NRTN_HUMAN
Immunogen	Recombinant human NRTN protein produced using CHO-based proprietary suspension cell line. Protein was purified from the cell culture supernatant. human Neurturin
Immunogen Region	
Specificity	Human NRTN
Immunogen Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081