

Anti-Neurturin antibody [3D12] (STJA0002882)

STJA0002882

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Neurturin is suitable for use in ELISA, Western Blot and Immunofluorescence research applications.
Applications	ELISA, WB, IF
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	3D12
Concentration	
Conjugation	Unconjugated
Purification	Protein G purification
Dilution Range	WB 0.2 to 1 µg/ml; ELISA 0.05 to 0.1 µg/ml; IF 0.5 to 5 µg/ml Monoclonal antibody working titer has to be established practically for each particular antigen and assay format
Formulation	PBS pH 7.4, with 0.1% sodium azide
Isotype	IgG1
Storage	Store at -20 or -70 °C upon receipt. Divide antibody into aliquots prior usage. Avoid multiple freeze-thaw cycles as product degradation may result
Instruction	

TARGET INFORMATION

Gene ID	4902
Gene Symbol	NRTN
Uniprot ID	NRTN_HUMAN
Immunogen	Recombinant human NRTN protein produced using CHO-based proprietary suspension cell line. Protein was purified from the cell culture supernatant. human Neurturin
Immunogen Region	
Specificity	Human NRTN
Immunogen Sequence	