

Anti-Recombinant-SARS-CoV-2 Nucleocapsid protein antibody [22F9] (STJA0002890)

STJA0002890

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Human monoclonal antibody anti-Recombinant-SARS-CoV-2 Nucleocapsid protein is suitable for use in ELISA and Western Blot research applications.
Applications	ELISA, WB
Host/Source	Human
Reactivity	SARS-CoV-2

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	22F9
Concentration	
Conjugation	Unconjugated
Purification	Protein A affinity chromatography
Dilution Range	ELISA 0.08-0.16 µg/ml
Formulation	PBS pH 7.4
Isotype	IgG1 lambda
Storage Instruction	Store at -20 °C to -70 °C long term. Avoid multiple freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Recombinant SARS-CoV-2 Nucleocapsid protein produced by CHO-based proprietary suspension cell line. Nucleocapsid protein SARS-CoV-2
Immunogen Region	
Specificity	SARS-CoV-2 Nucleocapsid protein
Immunogen Sequence	