

**Anti-Recombinant-SARS-CoV-2 Spike RBD antibody [28H3] (STJA0002891)**

STJA0002891

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Human monoclonal antibody anti-Recombinant-SARS-CoV-2 Spike RBD is suitable for use in ELISA research applications.
<b>Applications</b>	ELISA
<b>Host/Source</b>	Human
<b>Reactivity</b>	SARS-CoV-2

**PRODUCT PROPERTIES**

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	28H3
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Protein A affinity chromatography following desalting
<b>Dilution Range</b>	Indirect ELISA 0.02 – 0.04 ng/ml
<b>Formulation</b>	PBS pH 7.4
<b>Isotype</b>	IgG1 lambda
<b>Storage Instruction</b>	Store at -20 °C to -70 °C long term. Avoid multiple freeze-thaw cycles.

**TARGET INFORMATION**

<b>Gene ID</b>	
<b>Gene Symbol</b>	
<b>Uniprot ID</b>	
<b>Immunogen</b>	SARS-CoV-2 virus (PBMCs isolated from patients recovered from SARS-CoV-2 infection, antibodies) Spike RBD SARS-CoV-2
<b>Immunogen Region</b>	
<b>Specificity</b>	SARS-CoV-2 Spike RBD
<b>Immunogen Sequence</b>	

---

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081