

Anti-Ribonuclease 7 antibody [4F9] (STJA0002901)

STJA0002901

GENERAL INFORMATION



Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Ribonuclease 7 is suitable for use in ELISA, Western Blot, Immunohistochemistry and Immunofluorescence research applications.
Applications	ELISA, WB, IHC, IF
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	4F9
Concentration	
Conjugation	Unconjugated
Purification	Protein G purification
Dilution Range	ELISA 0.02 to 0.05 µg/ml; WB 0.5 to 1 µg/ml; IF 1 to 10 µg/ml; IHC 5 to 20 µg/ml Monoclonal antibody optimal dilution has to be established practically for each antigen and assay format
Formulation	PBS pH 7.4, with 0.1 % sodium azide
Isotype	IgG1
Storage	Store at -20 °C to -70 °C upon receipt. Divide antibody into aliquots prior usage. Avoid multiple freeze-thaw cycles as product degradation may result
Instruction	

TARGET INFORMATION

Gene ID	84659
Gene Symbol	RNASE7
Uniprot ID	RNAS7_HUMAN
Immunogen	Recombinant human RNase 7 protein produced by CHO-based proprietary suspension cell line. Human Ribonuclease 7
Immunogen Region	
Specificity	Human Ribonuclease 7
Immunogen Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081