

Anti-Cardiac Troponin I antibody [1B12] (STJA0002903)

STJA0002903

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Cardiac Troponin I is suitable for use in research applications.
Applications	
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	1B12
Concentration	
Conjugation	Unconjugated
Purification	
Dilution Range	ELISA, 0.05-0.1 Mu g/ml; WB, 0.2 to 0.5 Mu g/ml
Formulation	
Isotype	
Storage Instruction	Store at -20 to -70°C upon receipt. Divide antibody into aliquots prior usage. Avoid multiple freeze-thaw cycles

TARGET INFORMATION

Gene ID	7137
Gene Symbol	TNNI3
Uniprot ID	TNNI3_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	