

## Anti-S100-A9 protein antibody [6G10] (STJA0002914)

STJA0002914

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-S100-A9 protein is suitable for use in ELISA research applications.
<b>Applications</b>	ELISA
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	6G10
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Protein G purification
<b>Dilution Range</b>	ELISA 1-0.2 µg/ml
<b>Formulation</b>	PBS with 50% glycerol
<b>Isotype</b>	IgG1
<b>Storage Instruction</b>	Store at -20 to -70°C upon receipt. Divide antibody into aliquots prior usage. Avoid multiple freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	6280
<b>Gene Symbol</b>	S100A9
<b>Uniprot ID</b>	S10A9_HUMAN
<b>Immunogen</b>	. Human S100-A9
<b>Immunogen Region</b>	
<b>Specificity</b>	Human S100-A9
<b>Immunogen Sequence</b>	