

Anti-Recombinant-SARS-CoV-2 Spike S2 antibody [1H11] (STJA0002921)

STJA0002921

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Human monoclonal antibody anti-Recombinant-SARS-CoV-2 Spike S2 is suitable for use in ELISA research applications.
Applications	ELISA
Host/Source	Human
Reactivity	SARS-CoV-2

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	1H11
Concentration	
Conjugation	Unconjugated
Purification	Protein A affinity chromatography following desalting
Dilution Range	
Formulation	PBS pH 7.4
Isotype	IgG1 lambda
Storage Instruction	Store at -20 °C to -70 °C long term. Avoid multiple freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	SARS-CoV-2 virus (PBMCs isolated from patients recovered from SARS-CoV-2 infection, antibodies)
Immunogen Region	
Specificity	SARS-CoV-2 Spike S2
Immunogen Sequence	