

## Anti-Recombinant-HEL Human IgG1 antibody (STJA0002922)

STJA0002922

### GENERAL INFORMATION

|                          |                                                                                                             |
|--------------------------|-------------------------------------------------------------------------------------------------------------|
| <b>Product Type</b>      | Primary antibodies                                                                                          |
| <b>Short Description</b> | Cho Cells monoclonal antibody anti-Recombinant-HEL Human IgG1 is suitable for use in research applications. |
| <b>Applications</b>      |                                                                                                             |
| <b>Host/Source</b>       | CHO cells                                                                                                   |
| <b>Reactivity</b>        | Chicken                                                                                                     |

### PRODUCT PROPERTIES

|                            |                                                                      |
|----------------------------|----------------------------------------------------------------------|
| <b>Clonality</b>           | Monoclonal                                                           |
| <b>Clone ID</b>            |                                                                      |
| <b>Concentration</b>       |                                                                      |
| <b>Conjugation</b>         | Unconjugated                                                         |
| <b>Purification</b>        | Hitrap Mabelect Sure + Superdex 200 Increase                         |
| <b>Dilution Range</b>      |                                                                      |
| <b>Formulation</b>         | PBS pH 7.4                                                           |
| <b>Isotype</b>             | hIgG1                                                                |
| <b>Storage Instruction</b> | Store at -20...-70°C upon receipt. Avoid multiple freeze-thaw cycles |

### TARGET INFORMATION

|                           |                                 |
|---------------------------|---------------------------------|
| <b>Gene ID</b>            |                                 |
| <b>Gene Symbol</b>        |                                 |
| <b>Uniprot ID</b>         |                                 |
| <b>Immunogen</b>          | .                               |
| <b>Immunogen Region</b>   |                                 |
| <b>Specificity</b>        | Lysozyme from chicken egg white |
| <b>Immunogen Sequence</b> |                                 |