

## Anti-Recombinant-Procalcitonin antibody [6H5] (STJA0002944)

STJA0002944

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Cho Cells monoclonal antibody anti-Recombinant-Procalcitonin is suitable for use in ELISA research applications.
<b>Applications</b>	ELISA
<b>Host/Source</b>	CHO cells
<b>Reactivity</b>	Human

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	6H5
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Protein A affinity chromatography following desalting
<b>Dilution Range</b>	
<b>Formulation</b>	PBS pH 7.4
<b>Isotype</b>	IgG1
<b>Storage Instruction</b>	Store at -20 °C to -70 °C long term. Avoid multiple freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	
<b>Gene Symbol</b>	
<b>Uniprot ID</b>	
<b>Immunogen</b>	Produced by CHO-based cell line (expressed by QMCF Technology).
<b>Immunogen Region</b>	
<b>Specificity</b>	human PCT
<b>Immunogen Sequence</b>	