

Anti-Recombinant-Procalcitonin antibody [14E8] (STJA0002945)

STJA0002945

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Cho Cells monoclonal antibody anti-Recombinant-Procalcitonin is suitable for use in ELISA research applications.
Applications	ELISA
Host/Source	CHO cells
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	14E8
Concentration	
Conjugation	Unconjugated
Purification	Protein A affinity chromatography following desalting.
Dilution Range	
Formulation	PBS pH 7.4
Isotype	IgG1
Storage Instruction	Store at -20 °C to -70 °C long term. Avoid multiple freeze-thaw cycles.

TARGET INFORMATION

Gene ID	796
Gene Symbol	CALCA
Uniprot ID	CALC_HUMAN
Immunogen	Produced by CHO-based cell line (expressed by QMCF Technology).
Immunogen Region	
Specificity	human PCT
Immunogen Sequence	