

## Anti-Recombinant-HEL Mouse IgG2a antibody (STJA0002952)

STJA0002952

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Cho Cells monoclonal antibody anti-Recombinant-HEL Mouse IgG2a is suitable for use in research applications.
<b>Applications</b>	
<b>Host/Source</b>	CHO cells
<b>Reactivity</b>	Chicken

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Hitrap Mabselect Sure + Superdex 200 Increase
<b>Dilution Range</b>	
<b>Formulation</b>	PBS pH 7.4
<b>Isotype</b>	
<b>Storage Instruction</b>	Store at -20...-70°C upon receipt. Avoid multiple freeze-thaw cycles

### TARGET INFORMATION

<b>Gene ID</b>	
<b>Gene Symbol</b>	
<b>Uniprot ID</b>	
<b>Immunogen</b>	.
<b>Immunogen Region</b>	
<b>Specificity</b>	Lysozyme from chicken egg white
<b>Immunogen Sequence</b>	