

Anti-Cd80 antibody [16-10A1] {EV450} (STJA0003881)

STJA0003881

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Armenian Hamster monoclonal antibody anti-T-lymphocyte activation antigen CD80 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Armenian Hamster
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	16-10A1
Concentration	
Conjugation	EV450
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	12519
Gene Symbol	Cd80
Uniprot ID	CD80_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	