

Anti-CD8A antibody [OKT-8] {ER780} (STJA0003914)

STJA0003914

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-T-cell surface glycoprotein CD8 alpha chain is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	OKT-8
Concentration	
Conjugation	ER780
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG2ak
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	925
Gene Symbol	CD8A
Uniprot ID	CD8A_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen	
Sequence	