

Anti-CD45/Ptpnc antibody [30-F11] {EV450} (STJA0003929)

STJA0003929

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Receptor-type tyrosine-protein phosphatase C is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	30-F11
Concentration	
Conjugation	EV450
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG2bk
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	19264
Gene Symbol	Ptpnc
Uniprot ID	PTPRC_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen	
Sequence	