

Anti-CD45/PTPRC antibody [HI30] {ER780} (STJA0003931)

STJA0003931

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Primary antibodies |
| Short Description | Mouse monoclonal antibody anti-Receptor-type tyrosine-protein phosphatase C is suitable for use in Flow Cytometry research applications. |
| Applications | FC |
| Host/Source | Mouse |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|----------------------------|--|
| Clonality | Monoclonal |
| Clone ID | HI30 |
| Concentration | |
| Conjugation | ER780 |
| Purification | |
| Dilution Range | |
| Formulation | PBS with 0.05% Proclin300, 1% BSA |
| Isotype | IgG1k |
| Storage Instruction | Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product. |

TARGET INFORMATION

| | |
|-------------------------|-----------------------------|
| Gene ID | 5788 |
| Gene Symbol | PTPRC |
| Uniprot ID | PTPRC_HUMAN |
| Immunogen | |
| Immunogen Region | |
| Specificity | |
| Immunogen | |
| Sequence | |