

Anti-CEACAM8 antibody [G10F5] {FITC} (STJA0003945)

STJA0003945

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Carcinoembryonic antigen-related cell adhesion molecule 8 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	G10F5
Concentration	
Conjugation	FITC
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgMk
Storage	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.
Instruction	

TARGET INFORMATION

Gene ID	1088
Gene Symbol	CEACAM8
Uniprot ID	CEAM8_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	