

Anti-CD25/Il2ra antibody [PC-61.5.3] {ER780} (STJA0003947)

STJA0003947

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Interleukin-2 receptor subunit alpha is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	PC-61.5.3
Concentration	
Conjugation	ER780
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG1k
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	16184
Gene Symbol	Il2ra
Uniprot ID	IL2RA_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	