

## Anti-CD25/Il2ra antibody [PC-61.5.3] {ER780} (STJA0003948)

STJA0003948

### GENERAL INFORMATION

|                          |  |
|--------------------------|--|
| <b>Product Type</b>      | Primary antibodies   |
| <b>Short Description</b> | Rat monoclonal antibody anti-Interleukin-2 receptor subunit alpha is suitable for use in Flow Cytometry research applications. |
| <b>Applications</b>      | FC   |
| <b>Host/Source</b>       | Rat  |
| <b>Reactivity</b>        | Mouse  |

### PRODUCT PROPERTIES

|                            |  |
|----------------------------|--|
| <b>Clonality</b>           | Monoclonal   |
| <b>Clone ID</b>            | PC-61.5.3  |
| <b>Concentration</b>       |  |
| <b>Conjugation</b>         | ER780  |
| <b>Purification</b>        |  |
| <b>Dilution Range</b>      |  |
| <b>Formulation</b>         | PBS with 0.05% Proclin300, 1% BSA  |
| <b>Isotype</b>             | IgG1k  |
| <b>Storage Instruction</b> | Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product. |

### TARGET INFORMATION

|                           |                             |
|---------------------------|-----------------------------|
| <b>Gene ID</b>            | <a href="#">16184</a>       |
| <b>Gene Symbol</b>        | <a href="#">Il2ra</a>       |
| <b>Uniprot ID</b>         | <a href="#">IL2RA_MOUSE</a> |
| <b>Immunogen</b>          |                             |
| <b>Immunogen Region</b>   |                             |
| <b>Specificity</b>        |                             |
| <b>Immunogen Sequence</b> |                             |