

Anti-IL10 antibody [JES3-9D7] {FITC} (STJA0003950)

STJA0003950

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-Interleukin-10 is suitable for use in IConfocal MicroscopyC research applications.
Applications	ICFC
Host/Source	Rat
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	JES3-9D7
Concentration	
Conjugation	FITC
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG1k
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	3586
Gene Symbol	IL10
Uniprot ID	IL10_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	