

Anti-CD3 antibody [G4.18] {PE/TR} (STJA0003968)

STJA0003968

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD3 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	G4.18
Concentration	
Conjugation	PE/TR
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG3k
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	25710
Gene Symbol	Cd3d
Uniprot ID	CD3D_RAT
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	