

Anti-PROM1 antibody [W6B3C1] {APC} (STJA0003979)

STJA0003979

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Prominin-1 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	W6B3C1
Concentration	
Conjugation	APC
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG1k
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	8842
Gene Symbol	PROM1
Uniprot ID	PROM1_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	