

Anti-CD7 antibody [124-1D1] {FITC} (STJA0004014)

STJA0004014

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-T-cell antigen CD7 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	124-1D1
Concentration	
Conjugation	FITC
Purification	
Dilution Range	
Formulation	PBS with 0.05% Proclin300, 1% BSA
Isotype	IgG1k
Storage Instruction	Recommend storing between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	924
Gene Symbol	CD7
Uniprot ID	CD7_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	