

## Human PC1 (Procollagen 1) CLIA Kit (STJC0000143)

STJC0000143

### GENERAL INFORMATION


|                     |                                                                                                                                  |
|---------------------|----------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type</b> | Kits                                                                                                                             |
| <b>Description</b>  | This human PC1 kit is a highly sensitive in-vitro chemiluminescent immunoassay for the measurement of trace amounts of analytes. |
| <b>Applications</b> | CLIA                                                                                                                             |
| <b>Reactivity</b>   | Human                                                                                                                            |

### PRODUCT PROPERTIES

|                        |                                                                                                                                                                                             |
|------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Detection</b>       | Chemiluminescence                                                                                                                                                                           |
| <b>Sensitivity</b>     | 9.38pg/mL                                                                                                                                                                                   |
| <b>Detection Limit</b> | 15.63~1000pg/mL                                                                                                                                                                             |
| <b>Assay Time</b>      | 3.5h                                                                                                                                                                                        |
| <b>Storage</b>         | If unopened the kit may be stored at 2-8 °C for up to 1 month. If the kit will not be used within 1 month, store the components separately, according to the component table in the manual. |
| <b>Instruction</b>     |                                                                                                                                                                                             |
| <b>Note</b>            | FOR RESEARCH USE ONLY (RUO).                                                                                                                                                                |

### TARGET INFORMATION

|                     |                                           |
|---------------------|-------------------------------------------|
| <b>Sample Types</b> | Serum, plasma and other biological fluids |
| <b>Gene Symbol</b>  |                                           |
| <b>Uniprot ID</b>   |                                           |



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081