

Human LGALS2 (Galectin-2) CLIA Kit (STJC0000503)

STJC0000503

GENERAL INFORMATION


| | |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------|
| Product Type | Kits |
| Short Description | This human LGALS2 kit is a highly sensitive in-vitro chemiluminescent immunoassay for the measurement of trace amounts of analytes. |
| Applications | CLIA |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Detection | Chemiluminescence |
| Sensitivity | 37.5pg/mL |
| Detection Limit | 62.5~4000pg/mL |
| Assay Time | 3.5h |
| Storage | If unopened the kit may be stored at 2-8 °C for up to 1 month. If the kit will not be used within 1 month, store the components separately, according to the component table in the manual. |
| Instruction | |
| Note | FOR RESEARCH USE ONLY (RUO). |

TARGET INFORMATION

| | |
|---------------------|-------------------------------------------|
| Sample Types | Serum, plasma and other biological fluids |
| Gene Symbol | LGALS2 |
| Uniprot ID | LEG2_HUMAN |



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081