

## Human THBS2 (Thrombospondin-2) CLIA Kit (STJC0000720)

STJC0000720

### GENERAL INFORMATION

<b>Product Type</b>	Kits
<b>Description</b>	This human THBS2 kit is a highly sensitive in-vitro chemiluminescent immunoassay for the measurement of trace amounts of analytes.
<b>Applications</b>	CLIA
<b>Reactivity</b>	Human

### PRODUCT PROPERTIES

<b>Detection</b>	Chemiluminescence
<b>Sensitivity</b>	37.5pg/mL
<b>Detection Limit</b>	62.5~4000pg/mL
<b>Assay Time</b>	3.5h
<b>Storage Instruction</b>	If unopened the kit may be stored at 2-8 °C for up to 1 month. If the kit will not be used within 1 month, store the components separately, according to the component table in the manual.
<b>Note</b>	FOR RESEARCH USE ONLY (RUO).

### TARGET INFORMATION

<b>Sample Types</b>	Serum, plasma and other biological fluids
<b>Gene Symbol</b>	<a href="#">THBS2</a>
<b>Uniprot ID</b>	<a href="#">TSP2_HUMAN</a>