

Human IL33 (Interleukin-33) CLIA Kit (STJC0000974)

STJC0000974

GENERAL INFORMATION


| | |
|--------------------------|---|
| Product Type | Kits |
| Short Description | This human IL33 kit is a highly sensitive in-vitro chemiluminescent immunoassay for the measurement of trace amounts of analytes. |
| Applications | CLIA |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Detection | Chemiluminescence |
| Sensitivity | 4.69pg/mL |
| Detection Limit | 7.81~500pg/mL |
| Assay Time | 3.5h |
| Storage Instruction | If unopened the kit may be stored at 2-8 °C for up to 1 month. If the kit will not be used within 1 month, store the components separately, according to the component table in the manual. |
| Note | FOR RESEARCH USE ONLY (RUO). |

TARGET INFORMATION

| | |
|---------------------|---|
| Sample Types | Serum, plasma and other biological fluids |
| Gene Symbol | IL33 |
| Uniprot ID | IL33_HUMAN |



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081