

## Rat CASR (Calcium Sensing Receptor) CLIA Kit (STJC0001571)

STJC0001571

### GENERAL INFORMATION


|                          |   |
|--------------------------|---|
| <b>Product Type</b>      | Kits  |
| <b>Short Description</b> | This rat CASR kit is a highly sensitive in-vitro chemiluminescent immunoassay for the measurement of trace amounts of analytes. |
| <b>Applications</b>      | CLIA  |
| <b>Reactivity</b>        | Rat   |

### PRODUCT PROPERTIES

|                            |   |
|----------------------------|---|
| <b>Detection</b>           | Chemiluminescence   |
| <b>Sensitivity</b>         | 18.75pg/mL  |
| <b>Detection Limit</b>     | 31.25~2000pg/mL   |
| <b>Assay Time</b>          | 3.5h  |
| <b>Storage Instruction</b> | If unopened the kit may be stored at 2-8 °C for up to 1 month. If the kit will not be used within 1 month, store the components separately, according to the component table in the manual. |
| <b>Note</b>                | FOR RESEARCH USE ONLY (RUO).  |

### TARGET INFORMATION

|                     |   |
|---------------------|---|
| <b>Sample Types</b> | Serum, plasma and other biological fluids |
| <b>Gene Symbol</b>  | <a href="#">Casr</a>                      |
| <b>Uniprot ID</b>   | <a href="#">CASR_RAT</a>                  |



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081