

Rat IFN Alpha 1 (Interferon Alpha 1) CLIA Kit (STJC0001754)

STJC0001754

GENERAL INFORMATION


| | |
|--------------------------|--|
| Product Type | Kits |
| Short Description | This Rat IFN Alpha 1 kit is a highly sensitive in-vitro chemiluminescent immunoassay for the measurement of trace amounts of analytes. |
| Applications | CLIA |
| Reactivity | Rat |

PRODUCT PROPERTIES

| | |
|------------------------|------------------------------|
| Detection | Chemiluminescence |
| Sensitivity | 4.69pg/mL |
| Detection Limit | 7.81~500pg/mL |
| Assay Time | |
| Storage | |
| Instruction | |
| Note | FOR RESEARCH USE ONLY (RUO). |

TARGET INFORMATION

| | |
|---------------------|---------------------------------|
| Sample Types | Serum, Plasma, Cell supernatant |
| Gene Symbol | |
| Uniprot ID | |



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081