

Rat TDGF1 (Teratocarcinoma Derived Growth Factor 1) CLIA Kit (STJC0001901)

STJC0001901

GENERAL INFORMATION


| | |
|--------------------------|--|
| Product Type | Kits |
| Short Description | This rat TDGF1 kit is a highly sensitive in-vitro chemiluminescent immunoassay for the measurement of trace amounts of analytes. |
| Applications | CLIA |
| Reactivity | Rat |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Detection | Chemiluminescence |
| Sensitivity | 9.38pg/mL |
| Detection Limit | 15.63~1000pg/mL |
| Assay Time | 3.5h |
| Storage Instruction | If unopened the kit may be stored at 2-8 °C for up to 1 month. If the kit will not be used within 1 month, store the components separately, according to the component table in the manual. |
| Note | FOR RESEARCH USE ONLY (RUO). |

TARGET INFORMATION

| | |
|---------------------|---|
| Sample Types | Serum, plasma and other biological fluids |
| Gene Symbol | |
| Uniprot ID | |



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081