

Human AMH (AMH) Sandwich ELISA Kit (STJE0000005)

STJE0000005

GENERAL INFORMATION

Product Type	Kits
Short Description	This AMH Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Human

PRODUCT PROPERTIES

Detection	Colorimetric
Sensitivity	22.55 pg/ml
Detection Limit	156-10000pg/ml
Assay Time	4.0 hours
Storage	Store in the fridge for up to 1 month. After first, use store remaining components according to components table on insert.
Instruction Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	
Gene	AMH
Symbol	
Uniprot ID	MIS_HUMAN