

## Human IL10 (IL10) Sandwich ELISA Kit (STJE0000007)

STJE0000007

### GENERAL INFORMATION

<b>Product Type</b>	Kits
<b>Short Description</b>	This IL10 Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
<b>Applications</b>	ELISA
<b>Reactivity</b>	Human

### PRODUCT PROPERTIES

<b>Detection</b>	Colorimetric
<b>Sensitivity</b>	14.2 pg/mL
<b>Detection Limit</b>	31.2-2000 pg/ml
<b>Assay Time</b>	4.0 hours
<b>Storage Instruction</b>	Store in the fridge for up to 1 month. After first, use store remaining components according to components table on insert.
<b>Note</b>	FOR RESEARCH USE ONLY (RUO).

### TARGET INFORMATION

<b>Sample Types</b>	serum, plasma, cell culture supernates and other biological fluids
<b>Gene</b>	<a href="#">IL10</a>
<b>Symbol</b>	
<b>Uniprot ID</b>	<a href="#">IL10_HUMAN</a>