

Human CCL2 (CCL2) Sandwich ELISA Kit (STJE0000029)

STJE0000029

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Kits |
| Short Description | This CCL2 Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human cell culture supernatant, serum and plasma (EDTA, citrate, heparin). |
| Applications | ELISA |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|-------------------------|---|
| Detection | Colorimetric |
| Sensitivity | 9.12pg/mL |
| Detection Limit | 31.2-2000pg/ml |
| Assay Time | 4.0 hours |
| Storage | Store in the fridge for up to 1 month. After first, use store remaining components according to components table on insert. |
| Instruction Note | FOR RESEARCH USE ONLY (RUO). |

TARGET INFORMATION

| | |
|---------------------|----------------------------|
| Sample Types | |
| Gene | CCL2 |
| Symbol | |
| Uniprot ID | CCL2_HUMAN |