

Human ADIPOQ (ADIPOQ) Sandwich ELISA Kit (STJE0000033)

STJE0000033

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Kits |
| Short Description | This ADIPOQ Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human cell culture supernatant, serum and plasma (EDTA, citrate, heparin). |
| Applications | ELISA |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|-------------------------|---|
| Detection | Colorimetric |
| Sensitivity | 23.35 pg/mL |
| Detection Limit | 4000-62.5pg/ml |
| Assay Time | 4.0 hours |
| Storage | Store in the fridge for up to 1 month. After first, use store remaining components according to components table on insert. |
| Instruction Note | FOR RESEARCH USE ONLY (RUO). |

TARGET INFORMATION

| | |
|---------------------|-----------------------------|
| Sample Types | |
| Gene | ADIPOQ |
| Symbol | |
| Uniprot ID | ADIPO_HUMAN |