

Human CCL5 (CCL5) Sandwich ELISA Kit (STJE0000038)

STJE0000038

GENERAL INFORMATION

Product Type	Kits
Short Description	This CCL5 Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Human

PRODUCT PROPERTIES

Detection	Colorimetric
Sensitivity	11.24 pg/mL
Detection Limit	2000-31.2pg/ml
Assay Time	4.0 hours
Storage Instruction	Store in the fridge for up to 1 month. After first, use store remaining components according to components table on insert.
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	
Gene	CCL5
Symbol	
Uniprot ID	CCL5_HUMAN