

Human CXCL2 (CXCL2) Sandwich ELISA Kit (STJE0000079)

STJE0000079

GENERAL INFORMATION

Product Type	Kits
Short Description	This CXCL2 Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Human

PRODUCT PROPERTIES

Detection	Colorimetric
Sensitivity	2.54 pg/mL
Detection Limit	7.8-500 pg/ml
Assay Time	4.0 hours
Storage Instruction	Store in the fridge for up to 1 month. After first, use store remaining components according to components table on insert.
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	
Gene	CXCL2
Symbol	
Uniprot ID	CXCL2_HUMAN