

Human Albumin (Albumin) Sandwich ELISA Kit (STJE0000083)

STJE0000083

GENERAL INFORMATION

Product Type	Kits
Short Description	This Albumin Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	human

PRODUCT PROPERTIES

Detection	Colorimetric
Sensitivity	0.64ng/mL
Detection Limit	1.56-100ng/mL
Assay Time	4.0 hours
Storage	Store in the fridge for up to 1 month. After first, use store remaining components according to components table on insert.
Instruction Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	serum, plasma, cell culture supernates
Gene Symbol	
Uniprot ID	