

Human EPOR (EPOR) Sandwich ELISA Kit (STJE0000114)

STJE0000114

GENERAL INFORMATION

Product Type	Kits
Short Description	This EPOR Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Human

PRODUCT PROPERTIES

Detection	Colorimetric
Sensitivity	34.1pg/ml
Detection Limit	78-5000pg/ml
Assay Time	4.0 hours
Storage Instruction	Store in the fridge for up to 1 month. After first, use store remaining components according to components table on insert.
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	cell culture supernates, Cell Lysates, Tissue homogenates, serum and plasma.
Gene	EPOR
Symbol	
Uniprot ID	EPOR_HUMAN