

Human, Swine, Canine INS (INS) Sandwich ELISA Kit (STJE0000118)

STJE0000118

GENERAL INFORMATION

Product Type	Kits
Short Description	This INS Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human, swine, canine cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Human, Swine, Canine

PRODUCT PROPERTIES

Detection	Colorimetric
Sensitivity	6.18 pg/mL
Detection Limit	1000-15.6pg/ml
Assay Time	4.0 hours
Storage	Store in the fridge for up to 1 month. After first, use store remaining components according to components table on insert.
Instruction	
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	
Gene	INS
Symbol	
Uniprot ID	INS_HUMAN