

Human BAFF/TNFSF13B/BlyS/CD257 (BAFF/TNFSF13B/BlyS/CD257) Sandwich ELISA Kit (STJE0000136)

STJE0000136

GENERAL INFORMATION

Product Type	Kits
Short Description	This BAFF/TNFSF13B/BlyS/CD257 Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Human

PRODUCT PROPERTIES

Detection	Colorimetric
Sensitivity	24.53 pg/mL
Detection Limit	62.5-4000pg/ml
Assay Time	4.0 hours
Storage Instruction	Store in the fridge for up to 1 month. After first, use store remaining components according to components table on insert.
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	serum, plasma, cell Lysis and other biological fluids
Gene Symbol	
Uniprot ID	