

Human ANXA6 (ANXA6) Sandwich ELISA Kit (STJE0000265)

STJE0000265

GENERAL INFORMATION

Product Type	Kits
Short Description	This ANXA6 Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Human

PRODUCT PROPERTIES

Detection	Colorimetric
Sensitivity	0.069 ng/mL
Detection Limit	0.156-10ng/ml
Assay Time	4.0 hours
Storage Instruction	Store in the fridge for up to 1 month. After first, use store remaining components according to components table on insert.
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	serum, plasma, cell Lysis and other biological fluids
Gene	ANXA6
Symbol	
Uniprot ID	ANXA6_HUMAN