

Human BNP (BNP) Sandwich ELISA Kit (STJE0000345)

STJE0000345

GENERAL INFORMATION

Product Type	Kits
Description	This BNP Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Human

PRODUCT PROPERTIES

Detection	Colorimetric
Sensitivity	0.44 ng/mL
Detection Limit	0.78-50ng/ml
Assay Time	4.0 hours
Storage	Store in the fridge for up to 1 month. After first, use store remaining components according to components table on insert.
Instruction Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	serum, plasma, cell culture supernates and other biological fluids
Gene	
Symbol	
Uniprot ID	