

Human CAV1 (CAV1) Sandwich ELISA Kit (STJE0000400)

STJE0000400

GENERAL INFORMATION

Product Type	Kits
Short Description	This CAV1 Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Human

PRODUCT PROPERTIES

Detection	Colorimetric
Sensitivity	25.19 pg/mL
Detection Limit	5000-78.12pg/ml
Assay Time	
Storage Instruction	Store in the fridge for up to 1 month. After first, use store remaining components according to components table on insert.
Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	
Gene	CAV1
Symbol	
Uniprot ID	CAV1_HUMAN