

Human CD14 (CD14) Sandwich ELISA Kit (STJE0000418)

STJE0000418

GENERAL INFORMATION

Product Type	Kits
Description	This CD14 Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human cell culture supernatant, serum and plasma (EDTA, citrate, heparin).
Applications	ELISA
Reactivity	Human

PRODUCT PROPERTIES

Detection	Colorimetric
Sensitivity	0.13 ng/mL
Detection Limit	0.23-15ng/ml
Assay Time	4.0 hours
Storage	Store in the fridge for up to 1 month. After first, use store remaining components according to components table on insert.
Instruction Note	FOR RESEARCH USE ONLY (RUO).

TARGET INFORMATION

Sample Types	cell culture supernates, cell Lysis and other biological fluids
Gene	CD14
Symbol	
Uniprot ID	CD14_HUMAN