

Human COL1a1 (COL1a1) Sandwich ELISA Kit (STJE0000507)

STJE0000507

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Kits |
| Short Description | This COL1a1 Sandwich ELISA Kit is an in-vitro enzyme-linked immunosorbent assay for the measurement of samples in human cell culture supernatant, serum and plasma (EDTA, citrate, heparin). |
| Applications | ELISA |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|-------------------------|---|
| Detection | Colorimetric |
| Sensitivity | 12.2pg/ml |
| Detection Limit | 31.2-2000pg/ml |
| Assay Time | 4.0 hours |
| Storage | Store in the fridge for up to 1 month. After first, use store remaining components according to components table on insert. |
| Instruction Note | FOR RESEARCH USE ONLY (RUO). |

TARGET INFORMATION

| | |
|---------------------|--|
| Sample Types | serum, plasma and other biological fluids. |
| Gene | COL1A1 |
| Symbol | |
| Uniprot ID | CO1A1_HUMAN |